





Рис. 3. Схема создания рекомбинантной плазмиды pNZ::spike и проверка экспрессии целевого гена:

a – стратегия конструирования экспрессионной плазмиды pNZ::spike;

b – электрофореграмма с фрагментами рестрикции, использованными в приготовлении лигирующей смеси для последующей трансформации бактерий *B. subtilis* 168 (M – маркер молекулярной массы фрагментов ДНК (Thermo Scientific; кат. № SM0311));

в – результат идентификации белков с RBD в клетках штамма *L. lactis* pNZ::spike

(1 – лизат клеток без индукции, 2 – лизат клеток после индукции низином,

3 – положительный контроль RBD, M – маркер молекулярной массы белка (Bio-Rad; кат. № 161-0376))

Fig. 3. The scheme of recombinant plasmid pNZ::spike construction and verification of target gene expression:

a – the strategy of expression plasmid pNZ::spike engineering;

b – electropherogram of restriction fragments used in the ligation with the following transformation of *B. subtilis* 168 cells (M – DNA fragment molecular weight marker (Thermo Scientific; catalog No. SM0311));

c – result of identification of proteins with RBD in *L. lactis* pNZ::spike cells

(1 – cell lysate without induction, 2 – nisin-induced cell lysate,

3 – RBD positive control, M – protein molecular weight marker (Bio-Rad; catalog No. 161-0376))